# productinformation

## tesa® 62580

## Hydrophilic film for IVD products

tesa® 62580 is a transparent PET film functionalized on one side with a non-adhesive hydrophilic coating.

tesa® 62580 features:

- \*Fast and reliable fluid flow
- \*Proven compatibility with many enzymes
- \*Excellent aging stability
- \*Reliable bonding of the hydrophilic side to tesa spacer tapes
- \*Proven performance and integrity of multi-layers with tesa spacer tapes after aging

### Main Application

Securing fluid flow in capillary cell in-vitro diagnostic applications.

#### Technical Data

•	Contact angle water <	20	•	Color	transparent
•	Surface Tension >	60		Thickness of tape	100 μm
•	Backing material	PET film		Type of adhesive	functional

#### Additional Information

Storage and handling instructions:

Shelf life: 48 months from date of coating.

Avoid contact with water, material is sensitive against hunidity. Protect against mechanical impacts.